

Kallikrein-7 (Ile30) Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP03222
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human Kallikrein-7 (Cleaved-Ile30)
Mol wt	25kD
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	WB, IHC
Concentration	1 mg/ml
Full name	Kallikrein-7
Synonyms	Kallikrein-7 ;Cleaved-Ile30; Kallikrein-7; hK7; KLK7

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

N/A

Recommended Dilution

WB: 1: 500 - 1: 2000

IHC: 1: 50 - 1: 300

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com