

KLK1/2/3 (Cleaved-Ile25) Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP03300
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human KLK1/2/3 (Cleaved-Ile25)
Mol wt	28710
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	KLK1/2/3
Synonyms	KLK1/2/3 ;Cleaved-Ile25; Kallikrein-2; EC 3.4.21.35; Glandular kallikrein-1; hGK-1; Tissue kallikrein-2

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

N/A

Recommended Dilution

WB: 1: 1000 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year