

## LC3 A/B Rabbit Polyclonal Antibody

### Description

Product type	Primary Antibody
Code	BT-AP03381
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human LC3 A/B
Mol wt	13310
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	LC3 A/B
Synonyms	LC3 A/B

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

### Background

Probably involved in formation of autophagosomal vacuoles (autophagosomes).|PTM:The precursor molecule is cleaved by APG4B/ATG4B to form the cytosolic form| LC3-I. This is activated by APG7L/ATG7| transferred to ATG3 and conjugated to phospholipid to form the membrane-bound form| LC3-II.|Belongs to the MAP1 LC3 family.|subcellular location:LC3-II binds to the autophagic membranes.|subunit:3 different light chains| LC1| LC2 and LC3| can associate with MAP1A and MAP1B proteins.|tissue specificity:Most abundant in heart| brain| liver| skeletal muscle and testis but absent in thymus and peripheral blood leukocytes.|

### Recommended Dilution

WB: 1: 1000 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

### Images

No images.

### Storage

-20°C for 1 year