

## Galactosidase Beta1L3 Polyclonal Antibody

### Description

Product type	Primary Antibody
Code	BT-AP03446
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	The antiserum was produced against synthesized peptide derived from human BGLL3. AA range:291-340
Mol wt	69692
Species reactivity	Human
Clonality	Polyclonal
Recommended application	IHC-p, ELISA
Concentration	1 mg/ml
Full name	Galactosidase beta1L3 Antibody
Synonyms	GLB1L3; Beta-galactosidase-1-like protein 3

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

### Background

GLB1L3 (Galactosidase Beta 1 Like 3) is a Protein Coding gene. Diseases associated with GLB1L3 include leber congenital amaurosis. GO annotations related to this gene include hydrolase activity, hydrolyzing O-glycosyl compounds and beta-galactosidase activity. An important paralog of this gene is GLB1L2.


### Recommended Dilution

IHC: 1: 100 - 1: 300

ELISA: 1: 40000

Not yet tested in other applications.

### Images



Immunohistochemistry analysis of paraffin-embedded human colon carcinoma, using GLB1L3 Antibody. The picture on the right is blocked with the synthesized peptide.

### Storage

-20°C for one year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: [save@bt-laboratory.com](mailto:save@bt-laboratory.com) | [www.bt-laboratory.com](http://www.bt-laboratory.com)