

MPPD2 Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP03664
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human MPPD2
Mol wt	32340
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	Metallophosphoesterase MPPED2
Synonyms	Metallophosphoesterase MPPED2 ;EC 3.1.-.-;Fetal brain protein 239;239FB;Metallophosphoesterase domain-containing protein 2

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background


N/A

Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

Images



Western Blot analysis of Jurkat mouse-brain, rat-brain cell lysis. Primary Antibody was diluted at 1:1000. Secondary antibody was diluted at 1:10000

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com