

## GRTP1 Polyclonal Antibody

### Description

<b>Product type</b>	Primary Antibody
<b>Code</b>	BT-AP03840
<b>Host</b>	Rabbit
<b>Isotype</b>	IgG
<b>Size</b>	20ul, 50ul, 100ul
<b>Immunogen</b>	The antiserum was produced against synthesized peptide derived from human GRTP1. AA range:71-120
<b>Mol wt</b>	38557
<b>Species reactivity</b>	Human
<b>Clonality</b>	Polyclonal
<b>Recommended application</b>	WB, IHC-p, ELISA
<b>Concentration</b>	1 mg/ml
<b>Full name</b>	GRTP1 Antibody
<b>Synonyms</b>	GRTP1; TBC1D6; Growth hormone-regulated TBC protein 1; TBC1 domain family member 6

**This product is for research use only, not for use in human, therapeutic or diagnostic procedure.**

### Background

GRTP1 is a Protein Coding gene. Growth Hormone Regulated TBC Protein 1 may act as a GTPase-activating protein for Rab family protein(s). Gene Ontology (GO) annotations related to this gene include GTPase activator activity. An important paralog of this gene is TBC1D10C.

### Recommended Dilution


WB: 1: 500 - 1: 2000

IHC: 1: 100 - 1: 300


ELISA: 1: 40000

Not yet tested in other applications.


### Images



Western blot analysis of lysates from 293 cells, using GRTP1 Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using GRTP1 Polyclonal Antibody



Western blot analysis of the lysates from Jurkat cells using GRTP1 antibody.

### Storage

-20°C for one year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: [save@bt-laboratory.com](mailto:save@bt-laboratory.com) | [www.bt-laboratory.com](http://www.bt-laboratory.com)