

## NUD15 Rabbit Polyclonal Antibody

### Description

<b>Product type</b>	Primary Antibody
<b>Code</b>	BT-AP03875
<b>Host</b>	Rabbit
<b>Isotype</b>	IgG
<b>Size</b>	100ul, 50ul, 20ul
<b>Immunogen</b>	Synthesized peptide derived from human NUD15
<b>Mol wt</b>	18040
<b>Species reactivity</b>	Human, Mouse
<b>Clonality</b>	Polyclonal
<b>Recommended application</b>	WB
<b>Concentration</b>	1 mg/ml
<b>Full name</b>	Probable 8-oxo-dGTP diphosphatase NUDT15
<b>Synonyms</b>	Probable 8-oxo-dGTP diphosphatase NUDT15 ;8-oxo-dGTPase NUDT15;EC 3.6.1.55;7; 8-dihydro-8-oxoguanine-triphosphatase NUDT15;MutT homolog 2;MTH2;Nucleoside diphosphate-linked moiety X motif 15

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

### Background

N/A

### Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

### Images

No images.

### Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: [save@bt-laboratory.com](mailto:save@bt-laboratory.com) | [www.bt-laboratory.com](http://www.bt-laboratory.com)