

Numb(Phospho-Ser276) Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP03888
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized phospho peptide around human Numb (Ser276)
Mol wt	N/A
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	Numb
Synonyms	Numb ;Ser276 ; Protein numb homolog; h-Numb; Protein S171

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

The protein encoded by this gene plays a role in the determination of cell fates during development. The encoded protein whose degradation is induced in a proteasome-dependent manner by MDM2 is a membrane-bound protein that has been shown to associate with EPS15, LNX1 and NOTCH1. Alternative splicing results in multiple transcript variants.

Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com