

## OSB10 Rabbit Polyclonal Antibody

### Description

Product type	Primary Antibody
Code	BT-AP03946
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human OSB10
Mol wt	N/A
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	OSB10
Synonyms	OSB10; Oxysterol-binding protein-related protein 10; ORP-10; OSBP-related protein 10

**This product is for research use only, not for use in human, therapeutic or diagnostic procedure.**

### Background


This gene encodes a member of the oxysterol-binding protein (OSBP) family| a group of intracellular lipid receptors. Like most members| the encoded protein contains an N-terminal pleckstrin homology domain and a highly conserved C-terminal OSBP-like sterol-binding domain. Multiple transcript variants encoding different isoforms have been found for this gene.

### Recommended Dilution

WB: 1: 500 - 1: 2000

Not yet tested in other applications.

### Images



Western blot analysis of lysates from MDA-MB cells, primary antibody was diluted at 1:1000, 4°C overnight

### Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: [save@bt-laboratory.com](mailto:save@bt-laboratory.com) | [www.bt-laboratory.com](http://www.bt-laboratory.com)