

OSBP1 Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP03950
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human OSBP1
Mol wt	88770
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	OSBP1
Synonyms	OSBP1

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background


Oxysterol binding protein is an intracellular protein that is believed to transport sterols from lysosomes to the nucleus where the sterol down-regulates the genes for the LDL receptor| HMG-CoA reductase| and HMG synthetase

Recommended Dilution

WB: 1: 500 - 1: 2000

Not yet tested in other applications.

Images



Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com