

PHGDH(Phospho Thr204) Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP04133
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human PHGDH (phospho Thr204)
Mol wt	58630
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	PHGDH
Synonyms	PHGDH ;phospho Thr204; D-3-phosphoglycerate dehydrogenase; 3-PGDH; EC 1.1.1.95

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

N/A

Recommended Dilution

WB: 1: 1000 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com