

PRSS8 Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP04299
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human PRSS8 (heavy chain Cleaved-Ile45)
Mol wt	37730
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	PRSS8
Synonyms	PRSS8 ;heavy chain; Cleaved-Ile45; Prostatic; EC 3.4.21.-; Channel-activating protease 1; CAP1; Serine protease 8; Prostatic light chain; Prostatic heavy chain;

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

N/A

Recommended Dilution

WB: 1: 1000 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com