

PTGDS Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP04330
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human PTGDS
Mol wt	20900
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	Prostaglandin-H2 D-isomerase
Synonyms	Prostaglandin-H2 D-isomerase ;EC 5.3.99.2;Beta-trace protein;Cerebrin-28;Glutathione-independent PGD synthase;Lipocalin-type prostaglandin-D synthase;Prostaglandin-D2 synthase;PGD2 synthas

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhui Dist, Jiaying, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com