

RanBP3(Phospho-Ser58) Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP04384
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized phospho peptide around human RanBP3 (Ser58)
Mol wt	N/A
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	RanBP3
Synonyms	RanBP3 ;Ser58; Ran-binding protein 3; RanBP3

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

This gene encodes a protein with a RanBD1 domain that is found in both the nucleus and cytoplasm. This protein plays a role in nuclear export as part of a heteromeric complex. Alternate transcriptional splice variants| encoding different isoforms| have been characterized.

Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com