

RBM14 Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP04400
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human RBM14
Mol wt	73590
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	RNA-binding protein 14
Synonyms	RNA-binding protein 14 ;Paraspeckle protein 2;PSP2;RNA-binding motif protein 14;RRM-containing coactivator activator/modulator;Synaptotagmin-interacting protein;SYT-interacting protein

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com