

SLIB Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP04765
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human SLIB
Mol wt	11880
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	IHC-p, IF
Concentration	1 mg/ml
Full name	SLIB
Synonyms	SLIB

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

This gene encodes a member of the glucagon family of proteins. The encoded preproprotein is produced in the hypothalamus and cleaved to generate the mature factor| known as somatoliberin| which acts to stimulate growth hormone release from the pituitary gland. Variant receptors for somatoliberin have been found in several types of tumors| and antagonists of these receptors can inhibit the growth of the tumors. Defects in this gene are a cause of dwarfism| while hypersecretion of the encoded protein is a cause of gigantism. Alternative splicing results in multiple transcript variants| at least one of which encodes a preproprotein that is proteolytically processed.

Recommended Dilution

IHC-p: 1: 50 - 1: 200

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com