

KMT1B Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP04855
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	The antiserum was produced against synthesized peptide derived from human SUV39H2. AA range:111-160
Mol wt	46682
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	KMT1B Antibody
Synonyms	SUV39H2; KMT1B; Histone-lysine N-methyltransferase SUV39H2; Histone H3-K9 methyltransferase 2; H3-K9-HMTase 2; Lysine N-methyltransferase 1B; Suppressor of variegation 3-9 homolog 2; Su(var)3-9 homolo

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Histone-lysine N-methyltransferase SUV39H2 is a Histone methyltransferase that specifically trimethylates Lys-9 of histone H3 using monomethylated H3 Lys-9 as substrate. H3 Lys-9 trimethylation represents a specific tag for epigenetic transcriptional repression by recruiting HP1 (CBX1, CBX3 and/or CBX5) proteins to methylated histones. Mainly functions in heterochromatin regions, thereby playing a central role in the establishment of constitutive heterochromatin at pericentric and telomere regions. H3 Lys-9 trimethylation is also required to direct DNA methylation at pericentric repeats. SUV39H1 is targeted to histone H3 via its interaction with RB1 and is involved in many processes, such as cell cycle regulation, transcriptional repression and regulation of telomere length. May participate in regulation of higher-order chromatin organization during spermatogenesis. Recruited by the large PER complex to the E-box elements of the circadian target genes such as PER2 itself or PER1, contributes to the conversion of local chromatin to a heterochromatin-like repressive state through H3 Lys-9 trimethylation. Among its related pathways are Lysine degradation and Chromatin organization. Gene Ontology (GO) annotations related to this gene include chromatin binding and transcription regulatory region sequence-specific DNA binding. An important paralog of this gene is EHMT1.


Recommended Dilution

WB: 1: 500 - 1: 2000


ELISA: 1: 40000

Not yet tested in other applications.


Images



Western blot analysis of the lysates from 293 cells using SUV39H2 antibody.



Western Blot analysis of various cells using KMT1B Polyclonal Antibody cells nucleus.



Western Blot analysis of COLO205 cells using KMT1B Polyclonal Antibody cells nucleus.

Storage

-20°C for one year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com