

SQRD Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP04877
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Synthesized peptide derived from human SQRD
Mol wt	49500
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	SQRD
Synonyms	SQRD

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background


The protein encoded by this gene may function in mitochondria to catalyze the conversion of sulfide to persulfides| thereby decreasing toxic concentrations of sulfide. Alternative splicing results in multiple transcript variants that encode the same protein.

Recommended Dilution

WB: 1: 500 - 1: 2000

Not yet tested in other applications.

Images



Western blot analysis of lysates from AD293 cells, primary antibody was diluted at 1:1000, 4°C overnight

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com