

Lamin B1 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP04936
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Recombinant Protein of Lamin-B1
Mol wt	66408
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	Lamin B1 Antibody
Synonyms	LMNB1; LMN2; LMNB; Lamin-B1

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background


LMNB1 encodes one of the two B-type lamin proteins and is a component of the nuclear lamina. A duplication of LMNB1 is associated with autosomal dominant adult-onset leukodystrophy (ADLD). Alternative splicing results in multiple transcript variants.

Recommended Dilution

WB: 1: 1000 - 3000

Not yet tested in other applications.

Images



Western blot analysis of 1) 293T, 2) Mouse Brain, 3) Rat Brain tissue, diluted at 1:2000. Secondary antibody was diluted at 1:20000 cells nucleus.

Storage

-20°C for one year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com