

## TAF9B Rabbit Polyclonal Antibody

### Description

Product type	Primary Antibody
Code	BT-AP05016
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human TAF9B
Mol wt	27610
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	Transcription initiation factor TFIID subunit 9B
Synonyms	Transcription initiation factor TFIID subunit 9B ;Neuronal cell death-related protein 7;DN-7;Transcription initiation factor TFIID subunit 9-like;Transcription-associated factor TAFII31L

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

### Background

#### Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

### Images

No images.

### Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: [save@bt-laboratory.com](mailto:save@bt-laboratory.com) | [www.bt-laboratory.com](http://www.bt-laboratory.com)