

TDIF1 Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP05062
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Synthesized peptide derived from human TDIF1
Mol wt	36190
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	TDIF1
Synonyms	TDIF1

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background


DNTTIP1 binds DNA and enhances the activity of terminal deoxynucleotidyltransferase (TDT) or DNTT; MIM 187410| a DNA polymerase that catalyzes the polymerization of DNA in the absence of a DNA template (Yamashita et al. 2001)

Recommended Dilution

WB: 1: 500 - 1: 2000

Not yet tested in other applications.

Images



Western blot analysis of lysates from HEK293 cells, primary antibody was diluted at 1:1000, 4°C overnight

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com