

TPC1 Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP05221
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Synthesized peptide derived from human TPC1
Mol wt	89760
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	TPC1
Synonyms	TPC1

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background


Voltage-gated Ca(2+) and Na+ channels have 4 homologous domains| each containing 6 transmembrane segments| S1 to S6. TPCN1 is similar to these channels| but it has only 2 domains containing S1 to S6 (Ishibashi et al.| 2000

Recommended Dilution

WB: 1: 500 - 1: 2000

Not yet tested in other applications.

Images



Western blot analysis of lysates from HpeG2 cells, primary antibody was diluted at 1:1000, 4°C overnight

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com