

TRH Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP05244
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human TRH
Mol wt	26620
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	Prothyroliberin [Cleaved into: Thyroliberin
Synonyms	Prothyroliberin [Cleaved into: Thyroliberin ;Protirelin;TSH-releasing factor;Thyrotropin-releasing factor;TRF;Thyrotropin-releasing hormone;TRH]

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year