

MAZ Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP05254
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	The antiserum was produced against synthesized peptide derived from human MAZ. AA range:159-208
Mol wt	48608
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	MAZ Antibody
Synonyms	MAZ; ZNF801; Myc-associated zinc finger protein; MAZI; Pur-1; Purine-binding transcription factor; Transcription factor Zif87; ZF87; Zinc finger protein 801

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

MAZ (MYC Associated Zinc Finger Protein) is a Protein Coding gene. Diseases associated with MAZ include mental retardation, autosomal dominant 31. Among its related pathways are SIDS Susceptibility Pathways. GO annotations related to this gene include nucleic acid binding and chromatin binding. An important paralog of this gene is PATZ1.


Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 10000

Not yet tested in other applications.

Images



Western blot analysis of the lysates from HepG2 cells using MAZ antibody.

Storage

-20°C for one year