

TSC2(Phospho Ser1254) Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP05280
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human TSC2 (Phospho Ser1254)
Mol wt	198770
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	TSC2
Synonyms	TSC2 ;Phospho Ser1254; Tuberin; Tuberous sclerosis 2 protein

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

alternative products:Additional isoforms seem to exist. Experimental confirmation may be lacking for some isoforms|disease:Defects in TSC2 are a cause of lymphangioleiomyomatosis (LAM)

Recommended Dilution

WB: 1: 1000 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com