

TSLP Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP05285
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human TSLP
Mol wt	17490
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	TSLP
Synonyms	TSLP; Thymic stromal lymphopoietin

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Cytokine that induces the release of T cell-attracting chemokines from monocytes and| in particular| enhances the maturation of CD11c(+) dendritic cells. Can induce allergic inflammation by directly activating mast cells.|induction:Released by primary epithelial cells in response to certain microbial products| physical injury| or inflammatory cytokines.|subunit:Interacts with a receptor composed of CRLF2 and IL7R.|tissue specificity:Expressed in a number of tissues including heart| liver and prostate.|

Recommended Dilution

WB: 1: 1000 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com