

## Tubulin alpha Rabbit Polyclonal Antibody(Yeast)

### Description

|                                |  |
|--------------------------------|--|
| <b>Product type</b>            | Primary Antibody                                       |
| <b>Code</b>                    | BT-AP05321   |
| <b>Host</b>                    | Rabbit   |
| <b>Isotype</b>                 | IgG  |
| <b>Size</b>                    | 20ul, 50ul, 100ul                                      |
| <b>Immunogen</b>               | Purified recombinant Full length Tubulin alpha (Yeast) |
| <b>Mol wt</b>                  | N/A  |
| <b>Species reactivity</b>      | Yeast  |
| <b>Clonality</b>               | Polyclonal   |
| <b>Recommended application</b> | WB   |
| <b>Concentration</b>           | 1 mg/ml  |
| <b>Full name</b>               | Tubulin alpha  |
| <b>Synonyms</b>                | Tubulin alpha; Tubulin alpha                           |

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

### Background

N/A


### Recommended Dilution

WB: 1: 5000 - 1: 20000

ELISA: 1: 20000

Not yet tested in other applications.

### Images



Western Blot analysis of Yeast cells using Tubulin alpha rabbit pAb(Yeast). Secondary antibody was diluted at 1:20000

### Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: [save@bt-laboratory.com](mailto:save@bt-laboratory.com) | [www.bt-laboratory.com](http://www.bt-laboratory.com)