

MIA2 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP05404
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	The antiserum was produced against synthesized peptide derived from human MIA2. AA range:122-171
Mol wt	61684
Species reactivity	Human
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	MIA2 Antibody
Synonyms	MIA2; Melanoma inhibitory activity protein 2

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

May play a role in the pathophysiology of liver disease and may serve as a marker of liver damage.


Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 20000

Not yet tested in other applications.

Images



Western Blot analysis of HuvEc cells using MIA2 Polyclonal Antibody

Storage

-20°C for one year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com