

## MLF1 Polyclonal Antibody

### Description

Product type	Primary Antibody
Code	BT-AP05450
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Synthesized peptide derived from MLF1 . at AA range: 130-210
Mol wt	30627
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	MLF1 Antibody
Synonyms	MLF1; Myeloid leukemia factor 1; Myelodysplasia-myeloid leukemia factor 1

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

### Background

MLF1 encodes an oncoprotein which is thought to play a role in the phenotypic determination of hemopoetic cells. Translocations between MLF1 and nucleophosmin have been associated with myelodysplastic syndrome and acute myeloid leukemia. Multiple transcript variants encoding different isoforms have been found for this gene.


### Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 20000

Not yet tested in other applications.

### Images



Western blot analysis of MLF1 Antibody. The lane on the right is blocked with the MLF1 peptide.

### Storage

-20°C for one year