

mPR Gamma Polyclonal Antibody

Description

| | |
|-------------------------|---|
| Product type | Primary Antibody |
| Code | BT-AP05542 |
| Host | Rabbit |
| Isotype | IgG |
| Size | 20ul, 50ul, 100ul |
| Immunogen | The antiserum was produced against synthesized peptide derived from human MPRG. AA range:281-330 |
| Mol wt | 38014 |
| Species reactivity | Human, Mouse, Rat |
| Clonality | Polyclonal |
| Recommended application | WB, IF, ELISA |
| Concentration | 1 mg/ml |
| Full name | mPRgamma Antibody |
| Synonyms | PAQR5; MPRG; Membrane progesterin receptor gamma; mPR gamma; Progesterin and adipoQ receptor family member 5; Progesterin and adipoQ receptor family member V |

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

PAQR5 (Progesterin And AdipoQ Receptor Family Member 5) is a Protein Coding gene. GO annotations related to this gene include receptor activity and steroid binding. An important paralog of this gene is PAQR6.

Recommended Dilution


WB: 1: 500 - 1: 2000

IF: 1: 200 - 1: 1000

ELISA: 1: 5000

Not yet tested in other applications.

Images



Storage

-20°C for one year