

## mPR Gamma Polyclonal Antibody

### Description

Product type	Primary Antibody
Code	BT-AP05542
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	The antiserum was produced against synthesized peptide derived from human MPRG. AA range:281-330
Mol wt	38014
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	WB, IF, ELISA
Concentration	1 mg/ml
Full name	mPRgamma Antibody
Synonyms	PAQR5; MPRG; Membrane progesterin receptor gamma; mPR gamma; Progesterin and adipoQ receptor family member 5; Progesterin and adipoQ receptor family member V

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

### Background

PAQR5 (Progesterin And AdipoQ Receptor Family Member 5) is a Protein Coding gene. GO annotations related to this gene include receptor activity and steroid binding. An important paralog of this gene is PAQR6.

### Recommended Dilution


WB: 1: 500 - 1: 2000

IF: 1: 200 - 1: 1000

ELISA: 1: 5000

Not yet tested in other applications.

### Images



### Storage

-20°C for one year