

## ZN503 Rabbit Polyclonal Antibody

### Description

|                         |                                              |
|-------------------------|----------------------------------------------|
| Product type            | Primary Antibody                             |
| Code                    | BT-AP05696                                   |
| Host                    | Rabbit                                       |
| Isotype                 | IgG                                          |
| Size                    | 20ul, 50ul, 100ul                            |
| Immunogen               | Synthesized peptide derived from human ZN503 |
| Mol wt                  | 71060                                        |
| Species reactivity      | Human, Mouse                                 |
| Clonality               | Polyclonal                                   |
| Recommended application | WB                                           |
| Concentration           | 1 mg/ml                                      |
| Full name               | ZN503                                        |
| Synonyms                | ZN503                                        |

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

### Background


N/A

### Recommended Dilution

WB: 1: 500 - 1: 2000

Not yet tested in other applications.

### Images



Western blot analysis of lysates from PC-12 cells, primary antibody was diluted at 1:1000, 4°C overnight

### Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: [save@bt-laboratory.com](mailto:save@bt-laboratory.com) | [www.bt-laboratory.com](http://www.bt-laboratory.com)