

APM2 Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP05822
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human APM2
Mol wt	8360
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	IHC-p, IF
Concentration	1 mg/ml
Full name	APM2
Synonyms	APM2; C10orf116; Adipose specific 2

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

APM2 gene is exclusively expressed in adipose tissue. Its function is currently unknown.

Recommended Dilution

IHC-p: 1: 50 - 1: 200

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com