

Neurocan Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP05926
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Synthesized peptide derived from Neurocan . at AA range: 370-450
Mol wt	143093
Species reactivity	Human
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	Neurocan Antibody
Synonyms	NCAN; CSPG3; NEUR; Neurocan core protein; Chondroitin sulfate proteoglycan 3

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Neurocan is a chondroitin sulfate proteoglycan thought to be involved in the modulation of cell adhesion and migration.


Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 40000

Not yet tested in other applications.

Images



Western Blot analysis of Jk cells using Neurocan Polyclonal Antibody

Storage

-20°C for one year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com