

ADA23 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP06048
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human protein . at AA range: 720-800
Mol wt	N/A
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	Disintegrin and metalloproteinase domain-containing protein 23
Synonyms	Disintegrin and metalloproteinase domain-containing protein 23 ;ADAM 23;Metalloproteinase-like; disintegrin-like; and cysteine-rich protein 3;MDC-3

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

This gene encodes a member of the ADAM (a disintegrin and metalloprotease domain) family. Members of this family are membrane-anchored proteins structurally related to snake venom disintegrins and have been implicated in a variety of biological processes involving cell-cell and cell-matrix interactions| including fertilization| muscle development| and neurogenesis. It is reported that inactivation of this gene is associated with tumorigenesis in human cancers.

Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhui Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com