

Nkx-2.4 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP06051
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	The antiserum was produced against synthesized peptide derived from human NKX24. AA range:10-59
Mol wt	36011
Species reactivity	Human
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	Nkx-2.4 Antibody
Synonyms	NKX2-4; NKX2D; Homeobox protein Nkx-2.4; Homeobox protein NK-2 homolog D

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

NKX2-4 (NK2 Homeobox 4) is a Protein Coding gene. Gene Ontology (GO) annotations related to NKX2-4 include sequence-specific DNA binding. An important paralog of NKX2-4 is DLX5. UniProtKB/Swiss-Prot for NKX2-4 Gene.


Recommended Dilution

WB: 1: 500 - 1: 2000


ELISA: 1: 10000

Not yet tested in other applications.

Images



Western Blot analysis of various cells using Nkx-2.4 Polyclonal Antibody cells nucleus.



Western blot analysis of the lysates from K562 cells using NKX24 antibody.

Storage

-20°C for one year