

ASH1L Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP06090
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human ASH1L
Mol wt	326590
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	Histone-lysine N-methyltransferase ASH1L
Synonyms	Histone-lysine N-methyltransferase ASH1L; EC 2.1.1.43;ASH1-like protein;huASH1;Absent small and homeotic disks protein 1 homolog;Lysine N-methyltransferase 2H

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com