

ATG14L Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP06542
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Recombinant Protein of ATG14L
Mol wt	N/A
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	IHC-p, IF
Concentration	1 mg/ml
Full name	Beclin 1-associated autophagy-related key regulator
Synonyms	Beclin 1-associated autophagy-related key regulator ;Barkor;Autophagy-related protein 14-like protein;Atg14L; Beclin 1-associated autophagy-related key regulator; Barkor; Autophagy-related protein 14-like protein; Atg14L

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

miscellaneous:The sequence shown here is derived from an EMBL/GenBank/DDBJ third party annotation (TPA) entry.,

Recommended Dilution

IHC-p: 1: 50 - 1: 300

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com