

BMP4 Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP06580
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human BMP4
Mol wt	44880
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	BMP4
Synonyms	BMP4; Bone morphogenetic protein 4; BMP-4; Bone morphogenetic protein 2B; BMP-2B

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

disease:Defects in BMP4 are the cause of microphthalmia syndromic type 6 (MCOPS6)

Recommended Dilution

WB: 1: 1000 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com