

BRCA1(Phospho Ser988) Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP06585
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human BRCA1 (Phospho Ser988)
Mol wt	204930
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	BRCA1
Synonyms	BRCA1 ;Phospho Ser988; Breast cancer type 1 susceptibility protein; EC 6.3.2.-; RING finger protein 53

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

disease:Defects in BRCA1 are a cause of genetic susceptibility to breast cancer (BC)

Recommended Dilution

WB: 1: 1000 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com