

B3A2 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP06634
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human protein . at AA range: 830-910
Mol wt	N/A
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	Anion exchange protein 2
Synonyms	Anion exchange protein 2 ;AE 2;Anion exchanger 2;Non-erythroid band 3-like protein;BND3L;Solute carrier family 4 member 2

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

alternative products:Additional isoforms seem to exist,Plasma membrane anion exchange protein of wide distribution.,Belongs to the anion exchanger (TC 2.A.31) family.,

Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com