

c-Myb (Acetyl Lys480) Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP06664
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human c-Myb (Acetyl Lys480)
Mol wt	70400
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	c-Myb
Synonyms	c-Myb ; efg; Acetyl Lys480; Transcriptional activator Myb; Proto-oncogene c-Myb

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Comprised of 3 domains; an N-terminal DNA-binding domain| a centrally located transcriptional activation domain and a C-terminal domain involved in transcriptional repression.|Transcriptional activator; DNA-binding protein that specifically recognize the sequence 5'-YAAC

Recommended Dilution

WB: 1: 1000 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com