

## p130 Cas Polyclonal Antibody

### Description

Product type	Primary Antibody
Code	BT-AP06739
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	The antiserum was produced against synthesized peptide derived from human p130 Cas. AA range:131-180
Mol wt	93362
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	WB, IHC-p, ELISA
Concentration	1 mg/ml
Full name	p130 Cas Antibody
Synonyms	BCAR1; CAS; CASS1; CRKAS; Breast cancer anti-estrogen resistance protein 1; CRK-associated substrate; Cas scaffolding protein family member 1; p130cas

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

### Background

BCAR1 (BCAR1, Cas family scaffolding protein), or CAS, is an Src (MIM 190090) family kinase substrate involved in various cellular events, including migration, survival, transformation, and invasion.

### Recommended Dilution


WB: 1: 500 - 1: 2000

IHC: 1: 100 - 1: 300


ELISA: 1: 10000

Not yet tested in other applications.

### Images



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using p130 Cas Antibody. The picture on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using p130 Cas Polyclonal Antibody

### Storage

-20°C for one year