

CABP2 Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP06769
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human CABP2
Mol wt	24200
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	CABP2
Synonyms	CABP2; Calcium-binding protein 2; CaBP2

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

This gene belongs to a subfamily of calcium binding proteins that share similarity to calmodulin. Like calmodulin these family members can likely stimulate calmodulin-dependent kinase II and the protein phosphatase calcineurin. Calcium binding proteins are an important component of calcium mediated cellular signal transduction. Two transcript variants encoding different isoforms have been found for this gene.

Recommended Dilution

WB: 1: 1000 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com