

## p70 S6 kinase Alpha Polyclonal Antibody

### Description

Product type	Primary Antibody
Code	BT-AP06858
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Recombinant Protein of p70 S6 kinase $\alpha$
Mol wt	59140
Species reactivity	Mouse, Rat
Clonality	Polyclonal
Recommended application	IF, WB, IHC-p
Concentration	1 mg/ml
Full name	p70 S6 kinase alpha Antibody
Synonyms	RPS6KB1; STK14A; Ribosomal protein S6 kinase beta-1; S6K-beta-1; S6K1; 70 kDa ribosomal protein S6 kinase 1; P70S6K1; p70-S6K 1; Ribosomal protein S6 kinase I; Serine/threonine-protein kinase 14A; p70

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

### Background

RPS6KB1 (ribosomal protein S6 kinase B1) encodes a member of the ribosomal S6 kinase family of serine/threonine kinases. The encoded protein responds to mTOR (mammalian target of rapamycin) signaling to promote protein synthesis, cell growth, and cell proliferation. Activity of RPS6KB1 has been associated with human cancer. Alternatively spliced transcript variants have been observed. The use of alternative translation start sites results in isoforms with longer or shorter N-termini which may differ in their subcellular localizations. There are two pseudogenes for RPS6KB1 on chromosome 17.

### Recommended Dilution

WB: 1: 1000 - 2000

IF: 1: 50 - 200

Not yet tested in other applications.

### Images

No images.

### Storage

-20°C for one year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: [save@bt-laboratory.com](mailto:save@bt-laboratory.com) | [www.bt-laboratory.com](http://www.bt-laboratory.com)