

Pbx 4 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP06964
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Synthesized peptide derived from Pbx 4 . at AA range: 220-300
Mol wt	40854
Species reactivity	Human
Clonality	Polyclonal
Recommended application	IHC-p, ELISA
Concentration	1 mg/ml
Full name	Pbx 4 Antibody
Synonyms	PBX4; Pre-B-cell leukemia transcription factor 4; Homeobox protein PBX4

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

PBX4 (PBX homeobox 4) encodes a member of the pre-B cell leukemia transcription factor family. These proteins are homeobox proteins that play critical roles in embryonic development and cellular differentiation both as Hox cofactors and through Hox-independent pathways. The encoded protein contains a homeobox DNA-binding domain, but specific functions of the protein have not been determined. Alternatively spliced transcript variants have been observed for PBX4.

Recommended Dilution

IHC: 1: 100 - 1: 300

ELISA: 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for one year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com