

Peroxin 11Beta Polyclonal Antibody

Description

| | |
|--------------------------------|--|
| Product type | Primary Antibody |
| Code | BT-AP07070 |
| Host | Rabbit |
| Isotype | IgG |
| Size | 20ul, 50ul, 100ul |
| Immunogen | The antiserum was produced against synthesized peptide derived from human PEX11B. AA range:91-140 |
| Mol wt | 28431 |
| Species reactivity | Human, Mouse, Rat, Monkey |
| Clonality | Polyclonal |
| Recommended application | WB, ELISA |
| Concentration | 1 mg/ml |
| Full name | Peroxin 11beta Antibody |
| Synonyms | PEX11B; Peroxisomal membrane protein 11B; Peroxin-11B; Peroxisomal biogenesis factor 11B; Protein PEX11 homolog beta; PEX11-beta |

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Peroxisomal biogenesis factor 11 beta encoded by PEX11B facilitates peroxisomal proliferation and interacts with PEX19. Peroxisomal biogenesis factor 11 beta is found in the peroxisomal membrane. Several transcript variants, some protein-coding and some not protein-coding, have been found for PEX11B.


Recommended Dilution

WB: 1: 500 - 1: 2000


ELISA: 1: 10000

Not yet tested in other applications.

Images



Western blot analysis of lysates from COLO, HUVEC, COS7, and Jurkat cells, using PEX11B Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using Peroxin 11 β Polyclonal Antibody diluted at 1:500

Storage

-20°C for one year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com