

PHI-1 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP07120
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Synthesized peptide derived from PHI-1 . at AA range: 50-130
Mol wt	15911
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	IHC-p, ELISA
Concentration	1 mg/ml
Full name	PHI-1 Antibody
Synonyms	PPP1R14B; PLCB3N; PNG; Protein phosphatase 1 regulatory subunit 14B; Phospholipase C-beta-3 neighbouring gene protein

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Inhibitor of PPP1CA. Has over 50-fold higher inhibitory activity when phosphorylated.

Recommended Dilution

IHC: 1: 100 - 1: 300

ELISA: 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for one year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com