

Pim-2 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP07178
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Synthesized peptide derived from the C-terminal region of human Pim-2.
Mol wt	34190
Species reactivity	Human
Clonality	Polyclonal
Recommended application	IHC-p, ELISA
Concentration	1 mg/ml
Full name	Pim-2 Antibody
Synonyms	PIM2; Serine/threonine-protein kinase pim-2; Pim-2h

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

PIM2 encodes a protooncogene that acts as a serine/threonine protein kinase. Studies determined the encoded protein functions to prevent apoptosis and to promote cell survival.

Recommended Dilution

IHC: 1: 100 - 1: 300

ELISA: 1: 10000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for one year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com