

PPP1R2P3 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP07369
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Synthesized peptide derived from PPP1R2P3 . at AA range: 40-120
Mol wt	23048
Species reactivity	Human
Clonality	Polyclonal
Recommended application	IHC-p, ELISA
Concentration	1 mg/ml
Full name	PPP1R2P3 Antibody
Synonyms	PPP1R2P3; Putative protein phosphatase inhibitor 2-like protein 3; Protein phosphatase 1; regulatory subunit 2 pseudogene 3

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

PPP1R2P3 (Protein Phosphatase 1 Regulatory Inhibitor Subunit 2 Pseudogene 3) is a Pseudogene. GO annotations related to this gene include protein phosphatase inhibitor activity. nhibitor of protein-phosphatase 1.

Recommended Dilution

IHC: 1: 100 - 1: 300

ELISA: 1: 40000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for one year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com